

# PHLEBOTOMY & MULTI-USE CONTAINERS

# PATIENT ROOM SYSTEM

# WALL CABINET SYSTEM



**BEMIS PHLEBOTOMY AND MULTI-USE CONTAINERS** are designed for areas with restricted access.

DESCRIPTION	CATALOG NO.	PACKAGING
<b>One Quart Phlebotomy</b> Fits phlebotomy trays. Two separate openings.	100 030	100/case
<b>Phlebotomy Holder</b> Adhesive on bottom secures container to tabletop or tray.	400 000	20/case
<b>Two Gallon Container</b> Tabletop use or wall-mountable. Horizontal drop.	102 020 Beige 102 030 Red	30/case
<b>Three Gallon Container</b> Large capacity. Built-in handles for easy transport. Horizontal drop.	303 020 Beige 303 030 Red	12/case
<b>Mounting Bracket for Five Quart SharpSentinel® and Two Gallon Units</b> Sharps container slides on and locks into place.	435 020	10/case
<b>Key for Mounting Bracket</b> Unique plastic tool for removing units from mounting bracket #435.	440 020	10/case
<b>Three Gallon Wire Bracket (Non-locking)</b> Holds model 303 or 333 securely to wall.	423 000	5/case

**BEMIS PATIENT ROOM SYSTEM** is ideal for areas where space and staff supervision is limited.

DESCRIPTION	CATALOG NO.	PACKAGING
<b>Three Quart Midsize Wall Safe® and Five Quart Regular Wall Safe®</b> Rotating cylinder top allows one-handed operation. Locks onto concealed bracket.	125 020/150 020 Beige 125 030/150 030 Red	24/case
<b>Five Quart SharpSentinel®</b> Built-in handles for easy transport. Uses bracket #435, or use in wall cabinet.	175 020 Beige 175 030 Red	32/case
<b>Three Gallon SharpSentinel®</b> Large capacity. Wall-mounted or freestanding. Built-in handles.	333 020 Beige 333 030 Red	12/case
<b>Wall Safe Bracket</b> Securely locks Wall Safe in place. Concealed by container when in use.	410 020	10/case
<b>Wall Safe Key</b> Unique plastic tool for removal of midsize and regular Wall Safe units from bracket.	415 020	10/case
<b>Three Gallon Wire Bracket</b> Mounts easily to wall. Can be padlocked for added security (padlock not included).	425 000	5/case
<b>Padlock</b> For use with #425 locking bracket.	426 000	1/case

**BEMIS WALL CABINETS** secure sharps containers in a sturdy locked enclosure. Made of ABS plastic for durability and easy cleaning.

DESCRIPTION	CATALOG NO.	PACKAGING
<b>Five Quart Wall Cabinet</b> Fits sharps container model #175.	1051	1/case
<b>Three Gallon Wall Cabinet</b> Fits sharps container model #333.	1031	1/case
<b>Glove Box Holder</b> Snaps easily on to bottom of cabinet.	1041	1/case
<b>Shelf for Three Gallon Cabinet</b> Converts #1031 cabinet to fit #175 five quart sharps container.	1060	12/case

All Bemis sharps disposal containers meet or exceed OSHA, NIOSH and CDC requirements. Units are puncture resistant, exceeding current ASTM standards for wall thickness. All Bemis products are latex- and heavy metal-free for added safety.

# LARGE VOLUME & CHEMOTHERAPY CONTAINERS

# SHARPS DISPOSAL—CANADA

# SHARPSSENTINEL® SHARPS DISPOSAL SYSTEM



**BEMIS LARGE VOLUME AND CHEMOTHERAPY CONTAINERS** are ideal for high-volume areas and oncology units.

DESCRIPTION	CATALOG NO.	PACKAGING
Eight Gallon Container Built-in handles for easy transport. Dual openings.	108 030 Red	10/case
Eleven Gallon Container Accommodates large sharps and disposable instruments.	111 030 Red	6/case
Two Gallon Chemotherapy Container Screw-on cap. Gasketed lid for leak resistance.	202 004 White 202 040 Yellow	30/case
Eight Gallon Chemotherapy Container Built-in handles for easy transport. Dual openings. Gasketed lid.	208 004 White 208 040 Yellow	10/case
Eleven Gallon Chemotherapy Container Accepts larger amounts of chemotherapy waste. Gasketed lid.	211 004 White 211 040 Yellow	6/case
Rolling Cart Fits both eight- and eleven-gallon containers for convenient mobility.	430 030	1/case

**BEMIS SHARPS DISPOSAL UNITS FOR CANADA** feature appropriate yellow containers with English and French biohazard labels.

DESCRIPTION	CATALOG NO.	PACKAGING	DESCRIPTION	CATALOG NO.	PACKAGING
One Quart Container Fits phlebotomy trays. Dual openings.	100 040	100/case	Three Gallon Container For use in restricted access areas. Large opening lid.	303 040	12/case
Two Gallon Container For use in restricted access areas. Temporary and final closure.	102 040	30/case	Three Gallon SharpSentinel® Rotating cylinder top. Wall mounted or freestanding.	333 040	12/case
Wall Safe® Rotating lid allows one-handed disposal.	125 040 3 qt. 150 040 5 qt.	24/case 24/case	Eight and Eleven Gallon Containers For large volume disposal.	108 040 8 gal. 111 040 11 gal.	10/case 6/case
Five Quart SharpSentinel® Rotating cylinder top for hands-free operation.	175 040	32/case	Two Gallon Chemotherapy Container Screw-on cap. Gasketed lid for leak resistance.	202 040	30/case

See previous pages for hardware and bracket options.

